

Medical Gas Valve Box



We produce DFTT® Medical Gas Valve Box. We are a manufacturer. We devote to China medical machine with more time and energy and we are a good supplier. Our customers consist of Asian, Europe and so on. This Medical Gas Valve Box have powerful functions like containing almost all gas monitoring, using new technology and the best material. We would appreciate it if you can become our long term partner in your country.

Product Description

DFTT® **Medical Gas Valve Box** is made by Steel Plate, excellent quality, pipe joint thread of **Medical Gas Valve Box** is square shape, can be customized according to customer requirements. We sale **Medical Gas Valve Box** professionally, covering most countries around the world. We are expecting become your partner.

Product Parameter

Product name	Shape	material	size	color
Medical Gas valve box	Square Shape	Steel plate	52*62*15 cm	white customized

Medical Gas valve box Feature And Application





DFTT® **Medical Gas Valve Box** is used for Gas monitoring, often used in hospitals, where Medical staff control the Gas flow in each ward. **Medical Gas valve box** is a Medical Gas area valve box, including a box, a door, a lock and at least one valve set. Each valve set includes a sequentially connected inlet pipe, valve, intermediate pipe, barometer and outlet pipe. The front end of the box body is provided with a maintenance port, and the box door can be opened and closed with the maintenance port.

Medical Gas valve box Details



DFTT® **Medical Gas Valve Box** can customize the number of pipelines according to individual requirements, convenient and fast Gas inlet and outlet, and easy to adjust the Gas pressure, reduce the occurrence of accidents. **Medical Gas Valve Box** is easy to maintain, and has good after-sales service, a long warranty period, as well as a detailed explanation of the manual and more professional guidance. **Medical Gas Valve Box** built-in regulating valve, easy to adjust, suitable for all kinds of situations.